


Temprecord 5.28.0.2471 - [Logger S0061590]

File Program View Options Help

Logger S0061590 is waiting to be started. You may still alter the parameters if you wish.

 **Ready** to be started

```

Logger Type      : scientific
Serial Number    : S0061590
Last temperature : not available
Logging         : Temperature only
Sample Period    : 0:10:00 (6.0 / hour)
Start Delay      : 0:00:10
Start Time/Date  : OFF

Start with Button : ON
Stop with Button  : ON
Start with Switcher : OFF
Stop with Switcher : OFF

Allow Markers    : ON
Loop Overwrite   : ON
Safe Range Entry : OFF
Limit Delay      : 0

Upper Limit      : 28.00 °C
Lower Limit      : 18.00 °C


User Data       : Cherry's Sample Logger
:
:
:

Total Samples    : 138,456
Total Uses       : 5

Logger Version   : 2.14
T Sensor         : Standard
Memory           : 32k
  
```

Temprecord Options

- General
- COM Port
- File Saving
 - Web Files
 - Encryption
 - Date Formatting
- Export
 - Copy to Clipboard, Excel
- Summary
- Statistics
 - TTV
 - Refrigeration Index
 - Rate of Cooling
- Graph
 - Zoom and Presets
 - Values
 - View
- Parameter Defaults
 - Mon-T Defaults
 - Comment Labels
- Print/PDF
 - Print Colours
 - Headers and Footers
- Auto Mode
- Email
 - SMTP Email Options
- Web
 - Tracking
 - Proxy Settings
 - Web Load and Save
- Language
- System
- Diagnostics



Startup

Check for logger on startup

BlueTooth Compatibility

Lock toolbars

Units

Celsius

Fahrenheit

Prompts

<input checked="" type="checkbox"/> Prompt before starting	<input checked="" type="checkbox"/> Prompt before overwriting TR files
<input checked="" type="checkbox"/> Prompt before stopping	<input checked="" type="checkbox"/> Prompt before creating folders
<input checked="" type="checkbox"/> Prompt before reusing	<input checked="" type="checkbox"/> Prompt for PDF filename
<input checked="" type="checkbox"/> Prompt before exiting Temprecord	<input checked="" type="checkbox"/> Prompt before overwriting PDF
<input checked="" type="checkbox"/> Prompt on loggers with recalibration	<input checked="" type="checkbox"/> Prompt before overwriting export files
	<input checked="" type="checkbox"/> Prompt for Export Filename

OK

Cancel

Help